Sterilization in Personal Service Settings

| Sterilizer name: | Month: | Year: | |
|------------------|--------|-------|--|
| | | | |

Sterilization log

| Cycle date (day/month/year) | Load number | Contents | Start time | End time | Temperature | Pressure | Indicator change (yes or no) | Initials | Notes (maintenance, repairs, failed cycle, corrective actions) | Bi-weekly spore test (date and results) |
|--------------------------------|----------------|----------|---------------|-------------|-------------|----------|------------------------------------|----------|-----------------------------------------------------------------------|------------------------------------------------|
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |
| | | | | | | | | | | |

Please note:

- 1. Pressure cookers, glass-bead sterilizers, microwaves, U.V. sterilizers, boiling water, and dry-heat domestic ovens are not approved means of sterilizing.
- 2. The manufacturer's instructions for the sterilizing equipment must be onsite and easily accessible.
- 3. Records must be kept on-site for one year and available for an additional two years.
- 4. Failed spore tests must be reported to KFL&A Public Health.

